

### QUESTION 3: What methods can the Food and Drug Administration (FDA) and other regulatory bodies use to evaluate the efficacy of novel anti-infective technologies?

---

RECOMMENDATION: The FDA and other regulatory bodies can use *in vitro* cell culture methods to evaluate the antimicrobial efficacy against pathogens, followed by animal studies to evaluate osseointegration issues and a subsequent osteomyelitis/periprosthetic joint infection (PJI) animal model to evaluate the *in vivo* efficacy. However, clinical trials may be required for clearance or approval of some novel anti-infective technologies.

LEVEL OF EVIDENCE: Moderate

DELEGATE VOTE: Agree: 80%, Disagree: 3%, Abstain: 17% (Super Majority, Strong Consensus)

---

#### RATIONALE

Human clinical trials of anti-infective technologies are inherently difficult to perform according to Lazzarini et al. [1], due to the low incidence of implant-associated infections, the heterogeneous patient population, various treatment options in arthroplasty, the surrounding tissue condition after debridement and the broad range of causative pathogens and associated virulence patterns [2]. A cascade of *in vitro* cell culture methods and especially meaningful experimental animal models have to serve to fill this inevitable gap [1].

During the development of anti-infective biomaterials and devices and the determination of their anti-microbial properties, reliable *in vitro* test methods are essential to characterize implant surfaces [1,3]. In any evaluation procedure, cell proliferation has to be included as an important step in the course of infection [3]. For appropriate anti-microbial efficacy testing the independent aspects adhesion, proliferation and detection of bactericidal activity shall be considered in a consistent approach [3,4].

In the almost identical anti-microbial test methods, described with Japanese Industrial Standard (JIS) Z 2801:2010 and the International Organization for Standardization (ISO) 22196:2011 standards, the bacteria are applied onto the sample surface and covered under a sterile film, whereas for the American Society for Testing and Materials (ASTM) E 2180 test method the bacteria are applied as a thin agar slurry film. After 24 hours, by recovering vital bacteria from the samples, both test methods' anti-microbial efficacy is determined as the difference between the untreated reference and the anti-microbial sample. The major limitations are the required sample size (ISO 22196 5 x 5 cm, ASTM 3 x 3 cm) and the flat and smooth surface geometry, which is often not a given for orthopaedic implants [4]. In addition, hydrophobic surfaces can be unsuitable for testing according to ISO 22196, and the applied agar film (ASTM E 2180) can be too thick for non-leaching surface bound anti-microbials, thus leading to false-negative results.

Proliferation assay-based methods, first described by Bechert et al. [3], measure the antimicrobial efficacy based on the reproduction and release of daughter cells, monitoring the growth activity of these offspring bacteria over time. The main advantage of the proliferation-based assays is a broad applicability to flexible sample geometries (e.g., 2D and 3D), surface properties (e.g., smooth, textured, porous) and test conditions (e.g., leaching and non-leaching) [3–5]. Moreover, this method allows a parallelized investigation of many different setups in one test run ensuring a direct comparability, which results in increased explanatory power and higher sensitivity as given in the ISO and ASTM test methods [3,4]. However, the interpretation of test results is somehow more sophisticated, since growth of the offspring bacteria is analyzed rather than the vital cells on the sample surface [3,4]. In case of more complex surface structures and 3D geometries, which is the case for orthopaedic implants, the most reliable test method is a proliferation-based assay [4]. An important additional aspect is the contact of the implant to body fluids (such as blood, serum or interstitial liquid), having typically a high concentration of proteins, covering the device surface by a protein layer, which can have an impact on the antimicrobial performance of the material. Moreover, the influence of sterilization, aging degradation and persistence of the anti-microbial effect should be examined and testing should always be performed at least against gram-positive and gram-negative bacteria strains [4]. However, a direct transferability of *in vitro* results to *in vivo* performance is not stringently given. Thus, animal data are required to substantiate the antimicrobial efficacy *in vivo*.

To demonstrate unimpaired osseointegration for implant materials and surfaces that are modified by new anti-infective technologies in hip and knee arthroplasty, an appropriate animal study should be performed using controls based on long-term, clinically-established implant surfaces for cementless fixation, and also the base material and surface structure without the anti-infective treatment. Eto et al. [6] described a rat model with intramedullary implantation of a titanium rod to evaluate the osteoconductivity and osteogenesis in the meta- and diaphyseal region of the distal femur for experimental silver-oxide-containing hydroxyapatite coatings. They examined the implant anchorage strength at 2, 4 and 12 weeks post-implantation in a pull-out test, and performed a histological examination using a contralateral femur implantation with the same surface [6]. Analyzing the surface coverage with bone, they used this procedure to quantify the active peri-implant osteogenesis and osteoconductivity in the meta- and diaphysis of the femur in a comparison of anti-microbial surface treatments to a clinically-established hydroxyapatite (HA) coating [6]. Combining biomechanical and histological examinations, the model by Eto et al. [6] is valuable during the development phase of new anti-microbial implant surfaces to detect favorable solutions. The limitations of size, not allowing for testing multiple implants simultaneously and also significant dissimilarities between rat and human bone make a rat model unsuitable for clinically relevant osseointegration testing [7].

To evaluate new anti-microbial surface solutions for a clinical use in orthopaedic implants, their biocompatibility, peri-implant osteogenesis, osteoconductivity and ability of osseointegration should be tested in an animal model of a higher species, like sheep, goat, pig or dog [7,8]. Preferably a load-bearing model of the proximal tibia or distal femur in direct implantation site, or autologous left-right comparison should be performed, in reference to a clinically established surface (e.g., HA or porous coating) under a mid-term implantation duration of at least 26 weeks, to evaluate the osseointegration in a substantiated manner [7–10].

Animal models with osteomyelitis have been used previously to investigate potential treatment options using implants. After a review of the existing literature, it was found that a wide variety of osteomyelitis animal models exist [9]. However, no ideal single animal model exists to address implant associated osteomyelitis. Therefore, we propose that researchers and clinicians should ask indication and disease-specific questions and build on established appropriate animal models capable of answering their questions and enabling translations to the clinical situation [9]. Traditional methods to quantify bacterial load via colony forming unit (CFU) assays should be replaced with *in vivo* bio-luminescent imaging and radiological

outcome quantification. New anti-microbial treatments should be evaluated in regard to the host immune response utilizing biomarkers, and should be based on new technologies like the detection of bacteria by fluorescent in-situ hybridization in bone infection [9,11].

## REFERENCES

- [1] Lazzarini L, Overgaard KA, Conti E, Shirliff ME. Experimental osteomyelitis: what have we learned from animal studies about the systemic treatment of osteomyelitis? *J Chemother*. 2006;18:451–460. doi:10.1179/joc.2006.18.5.451.
- [2] Holinka J, Windhager R. [Management of prosthetic joint infections]. *Orthopade*. 2016;45:359–373; quiz 374. doi:10.1007/s00132–016–3247–8.
- [3] Bechert T, Steinrücke P, Guggenbichler JP. A new method for screening anti-infective biomaterials. *Nat Med*. 2000;6:1053–1056. doi:10.1038/79568.
- [4] Bruenke J, Roschke I, Agarwal S, Riemann T, Greiner A. Quantitative comparison of the antimicrobial efficiency of leaching versus nonleaching polymer materials. *Macromol Biosci*. 2016;16:647–654. doi:10.1002/mabi.201500266.
- [5] Kittinger C, Marth E, Windhager R, Weinberg AM, Zarfel G, Baumert R, et al. Antimicrobial activity of gentamicin palmitate against high concentrations of staphylococcus aureus. *J Mater Sci Mater Med*. 2011;22:1447–1453. doi:10.1007/s10856–011–4333–4.
- [6] Eto S, Miyamoto H, Shobuike T, Noda I, Akiyama T, Tsukamoto M, et al. Silver oxide-containing hydroxyapatite coating supports osteoblast function and enhances implant anchorage strength in rat femur. *J Orthop Res*. 2015;33:1391–1397. doi:10.1002/jor.22903.
- [7] Pearce AI, Richards RG, Milz S, Schneider E, Pearce SG. Animal models for implant biomaterial research in bone: a review. *Eur Cell Mater*. 2007;13:1–10.
- [8] Nuss KM, Auer JA, Boos A, von Rechenberg B. An animal model in sheep for biocompatibility testing of biomaterials in cancellous bones. *BMC Musculoskelet Disord*. 2006;7:67. doi:10.1186/1471–2474–7–67.
- [9] Reizner W, Hunter JG, O'Malley NT, Southgate RD, Schwarz EM, Kates SL. A systematic review of animal models for *Staphylococcus aureus* osteomyelitis. *Eur Cell Mater*. 2014;27:196–212.
- [10] Ignatius A, Peraus M, Schorlemmer S, Augat P, Burger W, Leyen S, et al. Osseointegration of alumina with a bioactive coating under load-bearing and unloaded conditions. *Biomaterials*. 2005;26:2325–2332. doi:10.1016/j.biomaterials.2004.07.029.
- [11] Alt V, Lips KS, Henkenbehrens C, Muhrer D, Oliveira Cavalcanti MC, Cavalcanti-Garcia M, et al. A new animal model for implant-related infected non-unions after intramedullary fixation of the tibia in rats with fluorescent in situ hybridization of bacteria in bone infection. *Bone*. 2011;48:1146–1153. doi:10.1016/j.bone.2011.01.018.

